

Product Datasheet

Claudin 5 Polyclonal Antibody FITC Conjugated (orb1541271)

Description	Claudin 5 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	O00501
Dilution Range	Flow-Cyt=1:50-200
Expiration Date	12 months from date of receipt.