

Product Datasheet

CD34 Polyclonal Antibody APC Conjugated (orb1541475)

Description	CD34 Polyclonal Antibody APC Conjugated
Species/Host	Rabbit
Conjugation	APC
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P28906
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.