

Product Datasheet

IL2RA/CD25 Polyclonal Antibody PE-Cy5 Conjugated (orb1541512)

Description	IL2RA/CD25 Polyclonal Antibody PE-Cy5 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy5
Tested Applications	FC
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P01590
Dilution Range	Flow-Cyt=1:50-200
Expiration Date	12 months from date of receipt.