

Product Datasheet

Cytomegalovirus pp65 Polyclonal Antibody Bio Conjugated (orb1541612)

Description	Cytomegalovirus pp65 Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	IF, IHC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q6SW59
Dilution Range	WB=1:50-200 IHC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.