

Product Datasheet

Beta-2-MG Mouse Monoclonal Antibody HRP Conjugated(B2E5) (orb1541743)

Description	Beta-2-MG Mouse Monoclonal Antibody HRP Conjugated(B2E5)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P61769
Dilution Range	Western blotting: 1:100-1000Immunohistochemistry1:100-500
Expiration Date	12 months from date of receipt.