

Product Datasheet

GPR55 Antibody FITC Conjugated (orb1541901)

Description	GPR55 Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9Y2T6
Dilution Range	FCM: 1:20-100 Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.