

Product Datasheet

ERK5 (Phospho-Ser731+Thr733) Antibody HRP Conjugated (orb1542117)

Description	ERK5 (Phospho-Ser731+Thr733) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q13164
Dilution Range	Western blotting: 1:100- 1000Immunohistochemistry1:100- 500
Expiration Date	12 months from date of receipt.