

Product Datasheet

GFAP (Phospho-Ser8) Polyclonal Antibody Cy5.5 Conjugated (orb1542127)

Description	GFAP (Phospho-Ser8) Polyclonal Antibody Cy5.5 Conjugated
Species/Host	Rabbit
Conjugation	Cy5.5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P14136
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.