

## Product Datasheet

# CHEK1(Phospho-Ser317) Antibody FITC Conjugated (orb1542133)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | CHEK1(Phospho-Ser317)<br>Antibody FITC Conjugated |
| <b>Species/Host</b>        | Rabbit  |
| <b>Conjugation</b>         | FITC  |
| <b>Tested Applications</b> | IF  |
| <b>Storage</b>             | Store at 4°C for 12 months.                       |
| <b>Note</b>                | For research use only                             |
| <b>Clonality</b>           | Recombinant                                       |
| <b>Source</b>              | Rabbit  |
| <b>Uniprot ID</b>          | <a href="#">O14757</a>                            |
| <b>Dilution Range</b>      | Immunofluorescence: 1:50-200                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |