

## Product Datasheet

### CARM1 Antibody HRP Conjugated (orb1542282)

<b>Description</b>	CARM1 Antibody HRP Conjugated
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC-P, WB
<b>Storage</b>	Store at 4°C for 12 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">Q86X55</a>
<b>Dilution Range</b>	Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500
<b>Expiration Date</b>	12 months from date of receipt.