

Product Datasheet

Acetyl-Histone H2B(K20) Antibody Cy5 Conjugated (orb1542334)

Description	Acetyl-Histone H2B(K20) Antibody Cy5 Conjugated
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P0C1H6
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.