

## Product Datasheet

# SGK1 (Phospho-Ser78) Antibody FITC Conjugated (orb1542394)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | SGK1 (Phospho-Ser78)<br>Antibody FITC Conjugated |
| <b>Species/Host</b>        | Rabbit   |
| <b>Conjugation</b>         | FITC   |
| <b>Tested Applications</b> | IF   |
| <b>Storage</b>             | Store at 4°C for 12 months.                      |
| <b>Note</b>                | For research use only                            |
| <b>Clonality</b>           | Recombinant                                      |
| <b>Source</b>              | Rabbit   |
| <b>Uniprot ID</b>          | <a href="#">O00141</a>                           |
| <b>Dilution Range</b>      | Immunofluorescence: 1:50-200                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |