

Product Datasheet

Rb P105 RB (Phospho-Ser807 + Ser811) Antibody FITC Conjugated (orb1542397)

Description	Rb P105 RB (Phospho-Ser807 + Ser811) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P06400
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.