

Product Datasheet

PLC gamma 2 (Phospho-Tyr759) Polyclonal Antibody Cy3 Conjugated (orb1542402)

Description	PLC gamma 2 (Phospho-Tyr759) Polyclonal Antibody Cy3 Conjugated
Species/Host	Rabbit
Conjugation	Cy3
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P16885
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.