

Product Datasheet

PERK(Phospho-Thr980) Antibody FITC Conjugated (orb1542409)

Description	PERK(Phospho-Thr980) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9Z2B5
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.