

Product Datasheet

Mst1(Phospho-Thr183) Mst2(Phospho-Thr180) Antibody FITC Conjugated (orb1542413)

Description	Mst1(Phospho-Thr183) Mst2(Phospho-Thr180) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q13043
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.