

Product Datasheet

IRF3 (Phospho-Ser396) Polyclonal Antibody Cy5 Conjugated (orb1542423)

Description	IRF3 (Phospho-Ser396) Polyclonal Antibody Cy5 Conjugated
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q14653
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.