

Product Datasheet

FOXO4 FOXO1 (Phospho-Ser197 Ser256) Antibody HRP Conjugated (orb1542431)

Description	FOXO4 FOXO1 (Phospho-Ser197 Ser256) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P98177
Dilution Range	Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500
Expiration Date	12 months from date of receipt.