

Product Datasheet

CSF1R (Phospho-Tyr723) Antibody HRP Conjugated (orb1542441)

Description	CSF1R (Phospho-Tyr723) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P07333
Dilution Range	Western blotting: 1:100- 1000Immunohistochemistry1:100- 500
Expiration Date	12 months from date of receipt.