

Product Datasheet

ERK1 2(Phospho-Thr202 + Tyr204) Antibody Cy3 Conjugated (orb1542450)

Description	ERK1 2(Phospho-Thr202 + Tyr204) Antibody Cy3 Conjugated
Species/Host	Rabbit
Conjugation	Cy3
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P27361
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.