

Product Datasheet

14-3-3 protein zeta delta (Phospho-Ser58) Antibody HRP Conjugated (orb1542454)

Description	14-3-3 protein zeta delta (Phospho-Ser58) Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Dilution Range	Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500
Expiration Date	12 months from date of receipt.