

## Product Datasheet

# JAK2(Phospho-Tyr1007+Tyr1008) Antibody HRP Conjugated (orb1542612)

<b>Description</b>	JAK2(Phospho-Tyr1007+Tyr1008) Antibody HRP Conjugated
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC-P, WB
<b>Storage</b>	Store at 4°C for 12 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">O60674</a>
<b>Dilution Range</b>	Western blotting: 1:100- 1000Immunohistochemistry1:100- 500
<b>Expiration Date</b>	12 months from date of receipt.