

## Product Datasheet

# Coxsackie Adenovirus Receptor Antibody FITC Conjugated (orb1542642)

<b>Description</b>	Coxsackie Adenovirus Receptor Antibody FITC Conjugated
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	IF
<b>Storage</b>	Store at 4°C for 12 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">P78310</a>
<b>Dilution Range</b>	FCM: 1:20-100 Immunofluorescence: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.