

Product Datasheet

CD13 APN ANPEN Antibody FITC Conjugated (orb1542966)

| | |
|----------------------------|--|
| Description | CD13 APN ANPEN Antibody FITC Conjugated |
| Species/Host | Rabbit |
| Conjugation | FITC |
| Tested Applications | IF |
| Storage | Store at 4°C for 12 months. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Rabbit |
| Uniprot ID | P15144 |
| Dilution Range | FCM: 1:20- 100Immunofluorescence: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |