

Product Datasheet

Glucocorticoid Receptor (Phospho-Ser226) Antibody FITC Conjugated (orb1543005)

Description	Glucocorticoid Receptor (Phospho-Ser226) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P04150
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.