

Product Datasheet

Giantin GOLGB1 Polyclonal Antibody Cy5 Conjugated (orb1543008)

Description	Giantin GOLGB1 Polyclonal Antibody Cy5 Conjugated
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q14789
Dilution Range	FCM:1:20-100IF:1:50-200
Expiration Date	12 months from date of receipt.