

Product Datasheet

Acetyl Coenzyme A Carboxylase alpha (Phospho-Ser80) Antibody FITC Conjugated (orb1543053)

Description	Acetyl Coenzyme A Carboxylase alpha (Phospho-Ser80) Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Dilution Range	Immunofluorescence: 1:50-200
Expiration Date	12 months from date of receipt.