

Product Datasheet

B7-H6 Polyclonal Antibody Cy5 Conjugated (orb1543074)

| | |
|----------------------------|---|
| Description | B7-H6 Polyclonal Antibody Cy5 Conjugated |
| Species/Host | Rabbit |
| Conjugation | Cy5 |
| Tested Applications | IF |
| Storage | Store at 4°C for 12 months. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Rabbit |
| Uniprot ID | Q68D85 |
| Dilution Range | IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |