

Product Datasheet

RAMP2 Antibody HRP Conjugated (orb1543150)

| | |
|----------------------------|---|
| Description | RAMP2 Antibody HRP Conjugated |
| Species/Host | Rabbit |
| Conjugation | HRP |
| Tested Applications | IHC-P, WB |
| Storage | Store at 4°C for 12 months. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Rabbit |
| Uniprot ID | O60895 |
| Dilution Range | Western blotting: 1:100-1000 Immunohistochemistry: 1:100-500 |
| Expiration Date | 12 months from date of receipt. |