

Product Datasheet

Adeno-Associated Virus Capsid Protein VP1 Antibody HRP Conjugated (orb1543413)

Description	Adeno-Associated Virus Capsid Protein VP1 Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	IHC-P, WB
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Dilution Range	Western blotting: 1:100- 1000Immunohistochemistry1:100- 500
Expiration Date	12 months from date of receipt.