

Product Datasheet

Cytomegalovirus pp65 Antibody FITC Conjugated (orb1543712)

Description	Cytomegalovirus pp65 Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Store at 4°C for 12 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q6SW59
Dilution Range	Immunofluorescence: 1:50- 200
Expiration Date	12 months from date of receipt.