

Product Datasheet

Anti-SARS-CoV-2 Nucleocapsid antibody(DM37), Rabbit mAb (orb1544105)

Description	Anti-SARS-CoV-2 Nucleocapsid antibody(DM37), Rabbit mAb
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.