

Product Datasheet

Anti-His Tag Antibody(DM1), Rabbit mAb (orb1544132)

Description	Anti-His Tag Antibody(DM1), Rabbit mAb
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Preservatives	0.1% Procline 300 Constituents: 50% Glycerol; PBS, pH 7.4; 0.1% rAlbumin
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.