

Product Datasheet

Anti-CD47 (magrolimab biosimilar) mAb (orb1544148)

| | |
|----------------------------|---|
| Description | Anti-CD47 (magrolimab biosimilar) mAb |
| Species/Host | Human |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |