


Product Datasheet

CD19 antibody (Biotin) (orb154467)

Description	Mouse monoclonal antibody conjugated to Biotin which recognizes CD19 that amplifies the signaling ca...
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, ICC
Immunogen	Human CCL (chronic lymphocytic leukemia) cells
Target	CD19
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P15391
Entrez	930
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface stainin...