


Product Datasheet

Cytokeratin 19 antibody (FITC) (orb154471)

Description	Mouse monoclonal antibody conjugated to FITC which recognizes Cytokeratin 19 that belongs to the int...
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	MCF-7 human breast adenocarcinoma cell line
Target	Cytokeratin 19
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P08727
Entrez	3880
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of MCF-7 cells stained using ...