

Product Datasheet

MAPK8/MAPK9/MAPK10 Rabbit Polyclonal Antibody (orb1545086)

Description	MAPK8/MAPK9/MAPK10 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P53779 , P45984 , P45983
Expiration Date	12 months from date of receipt.