


Product Datasheet

Human Gp130/IL6ST ELISA Kit (orb154573)

Description	Human Gp130/IL6ST ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L24-E619
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:200, concentration in serum and plasma is around 190-446 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	103537 MW
Uniprot ID	P40189
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Gp130/IL6ST PicoKine ELISA Kit sta...