


Product Datasheet

Human NT-4/NTF4/Neurotrophin-4 ELISA Kit (orb154574)

Description	Human NT-4/NTF4/Neurotrophin-4 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G81-A210
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	22427 MW
Uniprot ID	P34130
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human NT-4 PicoKine ELISA Kit standard C...