


Product Datasheet

Human CCL13/MCP4 ELISA Kit (orb154575)

Description	Human CCL13/MCP4 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-T98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	For serum, dilution ratio is 1:1 and concentration is around 40 pg/ml. For plasma, dilution ratio is 1:1 and concentration is around 100 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10986 MW
Uniprot ID	Q99616
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CCL13/MCP4 PicoKine ELISA Kit stan...