


Product Datasheet

Mouse CCL6/C10 ELISA Kit (orb154585)

Description	Mouse CCL6/C10 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G22-A116
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	12984 MW
Uniprot ID	P27784
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CCL6/C10
PicoKine
ELISA Kit
standa...