


## Product Datasheet

### Human TLR1/TIL ELISA Kit (orb154589)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human TLR1/TIL ELISA Kit  |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 156 pg/ml - 10,000 pg/ml  |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: S22-N578   |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml  |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).  |
| <b>MW</b>                  | 90291 MW  |
| <b>Uniprot ID</b>          | <b>Q15399</b>   |
| <b>Sensitivity</b>         | 10 pg/ml  |
| <b>Expiration Date</b>     | Please enquire.   |



Human TLR1  
PicoKine  
ELISA Kit  
standard c...