


Product Datasheet

Mouse TLR2/Toll-like receptor 2 ELISA Kit (orb154590)

Description	Mouse TLR2/Toll-like receptor 2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q25-L590
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	89449 MW
Uniprot ID	Q9QUN7
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
TLR2
PicoKine
ELISA Kit
standard
C...