


Product Datasheet

Mouse GDF-15 ELISA Kit (orb154591)

Description	Mouse GDF-15 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S189-A303
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 180 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	33225 MW
Uniprot ID	Q9Z0J7
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse GDF-15 PicoKine ELISA Kit standard...