

Product Datasheet

XL147 analogue(PI3K inhibitor X) (orb154629)

Description	XL147 analogue (SAR245408) is a selective and reversible class I PI3K inhibitor for wild type and mu...
Note	For research use only
Formula	$C_{21}H_{16}N_6O_2S_2$
MW	448.52
CAS Number	[956958-53-5]
SMILES	<chem>C(S(NC1C(NC2C=CC3C(C=2)=NSN=3)=NC2C(N=1)=CC=CC=2)(=O)=O)1=CC=C(C)C=C1</chem>
Expiration Date	12 months from date of receipt.