


## Product Datasheet

### PCI-32765 (Ibrutinib) (orb154698)

<b>Description</b>	PCI-32765 (Ibrutinib) is a potent and highly selective Btk inhibitor with IC50 of 0.5 nM.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>25</sub> H <sub>24</sub> N <sub>6</sub> O <sub>2</sub>
<b>MW</b>	440.5
<b>CAS Number</b>	[936563-96-1]
<b>SMILES</b>	<chem>C=CC(=O)N1CCC[C@H](C1)N2C3=C(C(=N2)C4=CC=C(C=C4)OC5=CC=CC=C5)C(=NC=N3)N</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of PCI-32765 (Ibrutin...)