

## Product Datasheet

# Recombinant Human Claudin 18.2 Protein, His Tag (orb1550292)

<b>Description</b>	Recombinant Human Claudin 18.2 Protein, His Tag
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/ug by LAL method.
<b>Target</b>	Claudin 18.2
<b>Note</b>	For research use only
<b>Purity</b>	> 95%by SEC-HPLC
<b>Source</b>	KMP2197, Recombinant Human Claudin 18.2 Protein is produced by Expi293 Cells expression system. The target protein is expressed with sequence (Met1-Val261) of Human Claudin 18.2 (Accession #P56856-2) fused with a His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.