

Product Datasheet

Recombinant Human LILRB1/CD85j Protein, His Tag (orb1550293)

Description	Recombinant Human LILRB1/CD85j Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	LILRB1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2196, Recombinant Human LILRB1/CD85j Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gly24-His Tag458) of Human LILRB1 (Accession #NP_006660.4) fused with an 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.