

## Product Datasheet

# Recombinant Mouse B7-H5/Gi24/VISTA Protein, His Tag (orb1550297)

<b>Description</b>	Recombinant Mouse B7-H5/Gi24/VISTA Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	B7H5
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP2192, Recombinant Mouse B7-H5/Gi24/VISTA Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Phe33 - Ala191) of Mouse VISTA (Accession #NP_083008.1) fused with an 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.