

Product Datasheet

Recombinant Mouse CTLA-4/CD152 Protein, His Tag (orb1550298)

Description	Recombinant Mouse CTLA-4/CD152 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CTLA4
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2191, Recombinant Mouse CTLA-4/CD152 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Glu36-Asp161) of Mouse CTLA4 (Accession #NP_033973.2) fused with an 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.